

XPRESSO ZERO Series



THE XPRESSO ZERO SERIES USES EVOLUTION'S PATENTED ROMAN-ARCH DESIGN FOR SUPERIOR STRENGTH AND STABILITY.

With the new padded seat cushion and wide back strap pad, users can sit in comfort and walk in style. The Xpresso Zero comes with an all-new braking system, allowing for effortless handle height adjustment and easy maintenance.

FEATURES

SEAT 1.5" thick padded seat cushion for extra comfort

BRAKE SYSTEM Concealed cable / Easy to adjust and maintain

COLLAPSIBLE BASKET Fully covered for increased privacy. Foldable & removable

WHEELS Lightweight 8" EVA tires / Offset rear forks for wider gait

WIDE BACK STRAP Extra wide for added comfort

HAND GRIPS Removable design for easy maintenance

COLORS Carbon Black, Shiraz Red, Sky Blue, Utopia Blue, Aquamarine Blue, Pearl White, Laurel Green, Coppery Brown, Fuchsia, Chrome Silver



SPECIFICATIONS

MINI

15 lbs*
17.5" Seat Width
18" Seat Height
25" Width x 23.5" Length
28.5" – 31.5" Handle Height
300 lbs Weight Capacity

REGULAR

15 lbs*
17.5" Seat Width
21" Seat Height
26" Width x 23.5" Length
31" – 35" Handle Height
300 lbs Weight Capacity

TALL

16 lbs*
17.5" Seat Width
24" Seat Height
27" Width x 24" Length
33" – 37" Handle Height
350 lbs Weight Capacity

OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Phone Holder
- Slow Down Brakes
- Eyeglass Holder

*Net Weight of walker without accessories, backstrap and folding basket